



List of performance reference compounds in EWL-PRC.

Phenanthrene - d10 – 2mg/EWL-PRC

Acenaphthene – d10 - 2mg/EWL-PRC

Chrysene-d12 - 2mg/EWL-PRC

Fluorene-d10 - 2mg/EWL-PRC

Benzo(e)pyrene-d12 - 2mg/EWL-PRC

C-13 Labeled

PCB1 – 50 ng/EWL-PRC

PCB8 - 50 ng/EWL-PRC

PCB37 - 50 ng/EWL-PRC

PCB54 - 50 ng/EWL-PRC

OCN – 200 ng/EWL-PRC

The concentrations are approximate! Should not be used for quantification.

Composition and concentrations of PRC mixture were changed in 2007 and it is valid only for the following Lot numbers:

160.407A

160.407B

160.407C

160.407D

160.407E